

Brillouin-Light-Scattering Spectroscopy

Benjamin Jungfleisch

SFB/TR49 Student Seminar
20th – 21th of July 2010



- Spin waves
- Brillouin Light Scattering (BLS)
 - Quantum mechanical picture
 - Conventional experimental setup
- Applications
 - Time-resolved BLS
 - Phase-resolved BLS
 - Wavevector-resolved BLS
- Summary

Spin waves

Landau-Lifshitz and Gilbert equation


masterformate durch Klicken bearbeiten

$$\frac{d\vec{M}}{dt} = -\gamma(\vec{M} \times \vec{H}_{eff}) + \frac{\alpha_G}{M_s} \left(\vec{M} \times \frac{d\vec{M}}{dt} \right)$$



Distinction between different energy contributions

- Exchange energy
(generated by twist of neighbored spins, short range interaction)
- Dipolar energy
(generated by magnetic poles in long-wavelength spin waves, long range interaction)



Brillouin-Light-Scattering

Raman-Scattering

- Optical phonons
- Molecule vibrations

Frequencies up to several THz


Grating spectrometer

Brillouin-Light-Scattering

- Acoustic phonons
- Spin waves

Frequencies below ~ 500 GHz

Tandem Fabry-Pérot interferometer



Annihilation (Anti-Stokes)

Stokes:

Anti-Stokes:

Energy conservation law:

$$\omega_S = \omega_i - \omega$$


$$\omega_S = \omega_i + \omega$$

Momentum conservation law:

$$k_S = k_i - k$$

$$k_S = k_i - k$$

Fabry-Pérot interferometer (FPI)



The functionality of a Fabry-Pérot interferometer is based on multi-beam interference on two plane-parallel surfaces.

The **phase difference**:
$$\Delta\varphi = \left(\frac{2\pi}{\lambda} \right) 2nl \cos(\theta)$$

Transmittance and reflectivity

The transmittance function is given by:

$$T = \frac{(1-R)^2}{1+R^2-2R\cos(\Delta\varphi)} = \frac{1}{1+F\sin^2(\Delta\varphi/2)},$$

where R is the reflectivity and $F = \frac{4R}{(1-R)^2}$ is the coefficient of the finesse F .

Maximum transmission occurs when

$$m = 1, 2, 3, \dots \qquad 2nl \cos(\theta) = m\lambda,$$


The free spectral range (FSR) is given by :

$$\Delta\lambda = \frac{\lambda_0^2}{2nl \cos \theta + \lambda_0} \approx \frac{\lambda_0^2}{2nl \cos \theta}$$

The finesse F is given by:


$$F = \frac{\Delta\lambda}{\delta\lambda} = \frac{\pi}{2 \arcsin(1/\sqrt{F})}$$

$$F \approx \frac{\pi \sqrt{F}}{2} = \frac{\pi \sqrt{R}}{1-R}$$



Our setup: $F \approx 110$

Fabry-Pérot interferometer



Maximal transmission: $l = m \frac{\lambda}{2}$

 optical bandpass

 frequency measurement


$$\lambda(\Delta l) = \lambda_R \left(1 + \frac{\Delta l}{l_R} \right)$$

Linear
relation:

$$\Delta \nu(\Delta l) = |\nu_R - \nu| = \left| \frac{c}{\lambda_R} \left(1 - \frac{1}{1 + \frac{\Delta l}{l_R}} \right) \right| = \left| \frac{c}{\lambda_R} \frac{\Delta l}{l_R + \Delta l} \right| \approx \left| \frac{c}{\lambda_R} \frac{\Delta l}{l_R} \right|$$


for IR= 0,3 mm: $\Delta \nu = 500$ GHz

Operation of the TFPI




Possibility to suppress higher orders.

Multi-pass tandem Fabry-Pérot interferometer (TFPI)



- High contrast: more than **1:10¹⁰**
- High spectral resolution in the **Sub-GHz-Regime (up to 50 MHz)**
- Accessible frequency range: **0,2 GHz – 500 GHz**

J. Sandercock, www.jrs-si.ch
H. Schultheiß, www.tfpdas.de



- optical resolution: up to 250 nm
- 2D scanning stage
- controlling sample position while measuring
- frequency range: 200 MHz – 1 THz
- spectral resolution: up to 50 MHz
- accuracy: **better than 20 nm**
- high reproducibility

TFPDAS 4
(www.tfpdas.de)
by H. Schultheiss

Applications

1) **Time resolution**


H. Schultheiss et al,
J. Phys. D **41**, 164017 (2008)

2) **Phase resolution**


A. A. Serga et al., APL **89**, 063506 (2006)
F. Fohr et al., RSI **80**, 043903 (2009)

3) **Wavevector resolution**


C. Sandweg et al,
Rev. Sci. Instrum. **81**, 073902 (2010)




2)



3)




Applications: Time-resolved BLS



Study the decay of spin waves
and the dissipation processes

H. Schultheiss et al,
J. Phys. D **41**, 164017 (2008)




Applications: Phase-resolved BLS

Inelastically scattered light contains phase information

But: BLS signal is proportional to laser intensity


▶ $I \propto x^2(t) = (A \cdot e^{i(\omega t + \phi)})^2 = A^2$ ▶ phase-information is lost!

Idea: interference!



running spin wave: wave has not reached the barrier yet

phase resolved scan:




$H_{\text{ext}} = 1846 \text{ Oe}$
 $\nu_{\text{RF}} = 7,125 \text{ GHz}$
 $t_{\text{puls}} = 200 \text{ ns}$

YIG-film (1,6 x 7,7 mm)

barrier (cut, 20 μm)

standard position scan:




standing spin wave: wave is reflected back and forms a standing wave

Applications: Wavevector-resolved BLS

No translational invariance in z direction


Only $k_{\parallel} = k \sin(\theta)$ conserved



Maximum wavevector: $k < 2.10^5 \text{ cm}^{-1}$

(Backscattering geometry: $\Delta k = 2 k_{\parallel}$)

Wavevector selection by rotating the sample




C. Sandweg et al,
 Rev. Sci. Instrum. 81, 073902 (2010)

Brillouin light scattering

study spin waves...

...wavevector
resolved

...time
resolved



...space
resolved

...frequency
resolved

...phase
resolved